

Newheek

Mobile X-ray Table With Bucky

Model : NK-T2(PG-N)

Website : www.xraybed.com

Email : service@newheek.com

Quick Details

Properties: Medical X-ray Equipments & Accessories

Brand Name: Newheek

Product Name: Simple Table With Bucky (Plexiglass Table)

Model Number: NK-T2(PG-N)

Place of Origin: Shandong, China (Mainland)

Table Surface: Plexiglass

Cassette Tray: Yes

Mobile: Yes

Table Surface Size: 2018mmx730mmx700mm

Cassette Tray Travel: 1000mm±10mm

Custom Made: Yes

Table Load: 150kg

Application: C arm, UC arm, Mobile x ray

Function: X ray radiology

Color: White

Product Application

Newheek medical x ray radiology table can be matched with different radiography equipment, such as c arm, mobile, uc-arm etc. It also can be used on vet x ray unit.

Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com

Product Outlook Pictures



Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com



Technical Parameter



Name	X Ray Radiology Table-with Cassette Tray
Model	NK-T2(PG-N)
Surface size	2018mmx730mmx700mm
Surface material	Plexiglass
Cassette tray travel	1000mm ± 10mm

Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com

Newheek Advantages:

Original manufacturer of x-ray machine accessories for more than 16 years.

- √ Customers could find all kinds of x-ray machine parts here.
- √ Offer on line technological support.
- √ Promise super product quality with best price and service.
- √ Support the third part inspection before delivery.
- √ Ensure the shortest delivery time.

Contact Us

Email : service@newheek.com

Hot Line: +86 536 8882360

FAQ

The common problems on the **X-ray Table**

1. How many types of the X-ray Table in NEWHEEK?

There are Medical X-ray Table and Veterinary X-ray Table.

2. What is the Payment Method ?

We could accept payment by TT, Western Union, Visa Credit Card, Paypal etc.

3. What is the surface material of the Table?

Plexiglass or Medium Density Fiberboard (MDF) .

4. How long is the Lead Time?

Within 15 workdays for x ray machine table after receipt of payment.

Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com