



Newheek



Six-Way Radiology X-ray Table

Model : NK-T4(CF-H)



Website : www.xraybed.com

Email : service@newheek.com



Quick Details

Properties: Medical X-ray Equipments & Accessories

Brand Name: Newheek

Product Name: 6-way floating radiology table

Model Number: NK-T4(CF-H)

Place of Origin: Shandong, China (Mainland)

Table Surface Material: Carbon Fiber

Cassette Tray: No

Mobile: Yes

Table Surface: 6 direction move

Table Surface Size: 2020mmx800mm

Surface Lock: Footswitch

Electrical: Yes

Custom Made: Yes

Table Load: 200kg

Application: C arm, UC arm, Mobile x ray

Function: X ray radiology

Color: White

Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com

Product Application

Newheek medical x ray radiology table can be matched with different radiography equipment, such as c arm, mobile, uc-arm etc. It also can be used on vet x ray unit.

Technical Parameter



Name	X Ray Radiology Table -Surface with 6 direction move
Model	NK-T4(CF-H)
Surface size	2020mmX800mm
Surface material	Carbon Fiber
Surface lock	foot switch
Surface movement	manual
Surface vertical travel	450mm±20mm
Surface horizontal travel	250mm±10mm
Table lifting travel	200mm (electric)
Table lowest position	705mm ± 10mm
Table highest position	905mm ± 10mm

Weifang Newheek Electronic Tech Co., Ltd.

Product Outlook Pictures



Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com



Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com

Newheek Advantages:

Original manufacturer of x-ray machine accessories for more than 16 years.

- √ Customers could find all kinds of x-ray machine parts here.
- √ Offer on line technological support.
- √ Promise super product quality with best price and service.
- √ Support the third part inspection before delivery.
- √ Ensure the shortest delivery time.

Contact Us

Email : service@newheek.com

Hot Line: +86 536 8882360

FAQ

The common problems on the **X-ray Table**

1. How many types of the X-ray Table in NEWHEEK?

There are Medical X-ray Table and Veterinary X-ray Table.

2. What is the Payment Method ?

We could accept payment by TT, Western Union, Visa Credit Card, Paypal etc.

3. What is the surface material of the Table?

Plexiglass or Carbon fiber.

4. How long is the Lead Time?

Within 15 workdays for x ray machine table after receipt of payment.

Weifang Newheek Electronic Tech Co., Ltd.

Add: No. 1 Factory Building of No.13426 of Yuqing East Street, Weifang, Shandong, China.
Website: www.xraybed.com Tel: 86 536 8882360 Email: service@newheek.com